

# Hill Hi-Lo Table

Since to 4450 €HT

is designed to serve the doctor who uses the ACTIVATOR® instrument or applies the "ACTIVATOR METHODS®" technique.\* Featuring a smooth transition between soft start-up and landing, the Hill Hi-Lo will not disturb the patient with any sudden stopping or bumping. The standing foot board automatically folds down as the table top is lowered. Patient size and comfort are taken into account with a wedge-shaped foot section that conforms to the angle of the leg (also available with flat cushion). The Hill Hi-Lo is available in 20 decorator colors; the base is available in beige, gray or black.

## Spécifications

Length: 6' • Width: 27"  
Optional Width: 30" (add 116€)  
Height: 34", optional 27"



\* (ACTIVATOR and ACTIVATOR METHODS are trademarks of Activator Methods, Inc. and are used for identification purposes only. Equipment made or distributed by Hill Laboratories is not made or approved by Activator Methods, Inc.)

# Forsco Labs

Exclusive Distributor - HILL Laboratories

149 avenue du Maine 75014 Paris, France Tel:  
+33 1 83 597 160 - Fax:+33 1 92 447 924  
contac@forscolabs.de • www.forscolabs.de

[WWW.FORSCOLABS.DE](http://WWW.FORSCOLABS.DE)